

SURGICAL POWER SINGLE-USE TOOLS



Single-use Perforator



No. of Blade	Head Diameter 4 mm	Head Diameter 6 mm	Head Diameter 9 mm	Head Diameter 12 mm
2	LT04C.DS	LT06C.DS	LT09C.DS	LT12C.DS
3	LTT09C.DS			

ATTACHMENT

1,500 rpm	ZSB2A	ZSB2A	ZSB2A	ZSB2A
1,500 rpm	ZSB3	ZSB3	ZSB3	ZSB3

Single-use Cutter

Tapered Straight




Tool Type	Head Diameter (mm)	Head Length 18 mm
	1.6	XD1816A.DS

ATTACHMENT	40,000 rpm	XJ6 XJH12F XJH12R
------------	------------	-------------------------

Single-use Round Carbide Bur

Small Bore: 2.38 mm




Tool Type	Diameter (mm)	Shaft Length 70 mm	Shaft Length 95 mm	Shaft Length 125 mm	Shaft Length 150 mm
	1.0	CR1006B70.DS	CR1006B95.DS	CR1006B125.DS	CR1006B150.DS
	2.0	CR2008B70.DS	CR2008B95.DS	CR2008B125.DS	CR2008B150.DS
	3.0	CR3010B70.DS	CR3010B95.DS	CR3010B125.DS	CR3010B150.DS
	4.0	CR4012B70.DS	CR4012B95.DS	CR4012B125.DS	CR4012B150.DS
	5.0	CR5012B70.DS	CR5012B95.DS	CR5012B125.DS	CR5012B150.DS
ATTACHMENT	80,000 rpm	MSB2A-W70	MSB2A-W95	MSB2A-W125	MSB2A-W150
	40,000 rpm	MSB1A-Z70	MSB1A-Z95	MSB1A-Z125	MSB1A-Z150
		MSB-FJ1A	MSB-FJ1B	MSB-FJ1C	MSB-FJ1D

Single-use Round Diamond Bur

Fine

Small Bore: 2.38 mm



Tool Type	Diameter (mm)	Shaft Length 70 mm	Shaft Length 95 mm	Shaft Length 125 mm	Shaft Length 150 mm
	1.0	DR1017B70.DS	DR1017B95.DS	DR1017B125.DS	DR1017B150.DS
	2.0	DR2017B70.DS	DR2017B95.DS	DR2017B125.DS	DR2017B150.DS
	3.0	DR3014B70.DS	DR3014B95.DS	DR3014B125.DS	DR3014B150.DS
	4.0	DR4012B70.DS	DR4012B95.DS	DR4012B125.DS	DR4012B150.DS
	5.0	DR5012B70.DS	DR5012B95.DS	DR5012B125.DS	DR5012B150.DS
	6.0	DR6012B70.DS	DR6012B95.DS	DR6012B125.DS	DR6012B150.DS
	7.0	DR7012B70.DS	DR7012B95.DS	DR7012B125.DS	DR7012B150.DS
ATTACHMENT	80,000 rpm	MSB2A-W70	MSB2A-W95	MSB2A-W125	MSB2A-W150
	40,000 rpm	MSB1A-Z70	MSB1A-Z95	MSB1A-Z125	MSB1A-Z150
		MSB-FJ1A	MSB-FJ1B	MSB-FJ1C	MSB-FJ1D

Single-use Round Stainless-steel Bur

Small Bore: 2.38 mm



Tool Type	Diameter (mm)	Shaft Length 70 mm	Shaft Length 95 mm	Shaft Length 125 mm	Shaft Length 150 mm
	2.0	SR2008B70.DS	SR2008B95.DS	SR2008B125.DS	SR2008B150.DS
	3.0	SR3010B70.DS	SR3010B95.DS	SR3010B125.DS	SR3010B150.DS
	4.0	SR4012B70.DS	SR4012B95.DS	SR4012B125.DS	SR4012B150.DS
	5.0	SR5012B70.DS	SR5012B95.DS	SR5012B125.DS	SR5012B150.DS
ATTACHMENT	80,000 rpm	MSB2A-W70	MSB2A-W95	MSB2A-W125	MSB2A-W150
	40,000 rpm	MSB1A-Z70	MSB1A-Z95	MSB1A-Z125	MSB1A-Z150
		MSB-FJ1A	MSB-FJ1B	MSB-FJ1C	MSB-FJ1D

Single-use Match-head Stainless-steel Bur

Small Bore: 2.38 mm



Tool Type	Diameter (mm)	Shaft Length 70 mm	Shaft Length 95 mm	Shaft Length 125 mm	Shaft Length 150 mm
	1.7	SM1702B70.DS	SM1702B95.DS	SM1702B125.DS	SM1702B150.DS
	2.2	SM2202B70.DS	SM2202B95.DS	SM2202B125.DS	SM2202B150.DS
	3.0	SM3002B70.DS	SM3002B95.DS	SM3002B125.DS	SM3002B150.DS
ATTACHMENT	80,000 rpm	MSB2A-W70	MSB2A-W95	MSB2A-W125	MSB2A-W150
	40,000 rpm	MSB1A-Z70	MSB1A-Z95	MSB1A-Z125	MSB1A-Z150
		MSB-FJ1A	MSB-FJ1B	MSB-FJ1C	MSB-FJ1D

Single-use Fluted Acorn Stainless-steel Bur


Small Bore: 2.38 mm



Tool Type	Diameter (mm)	Shaft Length 70 mm	Shaft Length 95 mm	Shaft Length 125 mm	Shaft Length 150 mm
	4.0	SOA4008B70.DS	SOA4008B95.DS	SOA4008B125.DS	SOA4008B150.DS
	6.0	SOA6008B70.DS	SOA6008B95.DS	SOA6008B125.DS	SOA6008B150.DS
ATTACHMENT	80,000 rpm	MSB2A-W70	MSB2A-W95	MSB2A-W125	MSB2A-W150
	40,000 rpm	MSB1A-Z70	MSB1A-Z95	MSB1A-Z125	MSB1A-Z150
		MSB-FJ1A	MSB-FJ1B	MSB-FJ1C	MSB-FJ1D

Small Bore: 2.38 mm



Tool Type	Diameter (mm)	Shaft Length 70 mm	Shaft Length 95 mm	Shaft Length 125 mm	Shaft Length 150 mm
	5.0	SSA5008B70.DS	SSA5008B95.DS	SSA5008B125.DS	SSA5008B150.DS
	7.0	SSA7008B70.DS	SSA7008B95.DS	SSA7008B125.DS	SSA7008B150.DS
ATTACHMENT	80,000 rpm	MSB2A-W70	MSB2A-W95	MSB2A-W125	MSB2A-W150
	40,000 rpm	MSB1A-Z70	MSB1A-Z95	MSB1A-Z125	MSB1A-Z150
		MSB-FJ1A	MSB-FJ1B	MSB-FJ1C	MSB-FJ1D

Single-use Round Diamond Micro Bur

Small Bore: 2.38 mm



Tool Type	Diameter (mm)	Angle (°)	Effective Length 90 mm	Effective Length 98 mm	Effective Length 115 mm
	1.0	0	DRM1017X90.DS	DRM1017X98.DS	DRM1017X115.DS
	2.0	0	DRM2017X90.DS	DRM2017X98.DS	DRM2017X115.DS
	1.0	20	DRM1017X9020.DS	DRM1017X9820.DS	DRM1017X11520.DS
	2.0	20	DRM2017X9020.DS	DRM2017X9820.DS	DRM2017X11520.DS
ATTACHMENT	80,000 rpm 40,000 rpm		XMSB2-W XMSB1-Z	XMSB2-W XMSB1-Z	XMSB2-W XMSB1-Z

Single-use Round Stainless-steel Micro Bur

Small Bore: 2.38 mm



Tool Type	Diameter (mm)	Angle (°)	Effective Length 90 mm	Effective Length 98 mm	Effective Length 115 mm
	1.0	0	SRM1006X90.DS	SRM1006X98.DS	SRM1006X115.DS
	2.0	0	SRM2008X90.DS	SRM2008X98.DS	SRM2008X115.DS
	1.0	20	SRM1006X9020.DS	SRM1006X9820.DS	SRM1006X11520.DS
	2.0	20	SRM2008X9020.DS	SRM2008X9820.DS	SRM2008X11520.DS
ATTACHMENT	80,000 rpm 40,000 rpm		XMSB2-W XMSB1-Z	XMSB2-W XMSB1-Z	XMSB2-W XMSB1-Z

Single-use Round Diamond Slim Bur

Slim Bore: 1.19 mm



Tool Type	Diameter (mm)	Shaft Length 100 mm	Shaft Length 105 mm	Shaft Length 111 mm	Shaft Length 155 mm
	1.0			DR1017G111.DS	
	2.0		DR2017G105.DS	DR2017G111.DS	
	3.0	DR3014G100.DS		DR3014G111.DS	DR3014G155.DS
	4.0	DR4012G100.DS		DR4012G111.DS	DR4012G155.DS
ATTACHMENT	80,000 rpm	MSB2-WX MSBF4580 MSBF451280	MSB2-WX MSBF4580 MSBF451280	MSB2-WX MSBF4580 MSBF451280	MSB2-WX MSBF45135 MSBF4512135 MSBF4510C135 MSBF4512C135

Slim Bore: 1.19 mm



Tool Type	Diameter (mm)	Shaft Length 160 mm	Shaft Length 175 mm	Shaft Length 180 mm	Shaft Length 195 mm
	2.0	DR2017G160.DS		DR2017G180.DS	
	3.0		DR3014G175.DS		DR3014G195.DS
	4.0		DR4012G175.DS		DR4012G195.DS
ATTACHMENT	80,000 rpm	MSB2-WX MSBF45135 MSBF4512135 MSBF4510C135 MSBF4512C135	MSB2-WX MSBF45155 MSBF45H155 MSBF4512155	MSB2-WX MSBF45155 MSBF45H155 MSBF4512155	MSB2-WX MSBF45175 MSBF45H175 MSBF4512175 MSBF4512C175

Single-use Round Stainless-steel Slim Bur

Slim Bore: 1.19 mm



Tool Type	Diameter (mm)	Shaft Length 100 mm	Shaft Length 105 mm	Shaft Length 111 mm	Shaft Length 155 mm
	1.0			SR1006G111.DS	
	2.0		SR2006G105.DS	SR2006G111.DS	
	3.0	SR3010G100.DS		SR3010G111.DS	SR3010G155.DS
	4.0	SR4012G100.DS		SR4012G111.DS	SR4012G155.DS
ATTACHMENT	80,000 rpm	MSB2-WX MSBF4580 MSBF451280	MSB2-WX MSBF4580 MSBF451280	MSB2-WX MSBF4580 MSBF451280	MSB2-WX MSBF45135 MSBF4512135 MSBF4510C135 MSBF4512C135

Slim Bore: 1.19 mm



Tool Type	Diameter (mm)	Shaft Length 160 mm	Shaft Length 175 mm	Shaft Length 180 mm	Shaft Length 195 mm
	2.0	SR2006G160.DS		SR2006G180.DS	
	3.0		SR3010G175.DS		SR3010G195.DS
	4.0		SR4012G175.DS		SR4012G195.DS
ATTACHMENT	80,000 rpm	MSB2-WX MSBF45135 MSBF4512135 MSBF4510C135 MSBF4512C135	MSB2-WX MSBF45155 MSBF45H155 MSBF4512155	MSB2-WX MSBF45155 MSBF45H155 MSBF4512155	MSB2-WX MSBF45175 MSBF45H175 MSBF4512175 MSBF4512C175

Single-use Match-head Diamond Slim Bur

Slim Bore: 1.19 mm



Tool Type	Diameter (mm)	Shaft Length 100 mm	Shaft Length 105 mm	Shaft Length 155 mm	Shaft Length 160 mm
	2.0		DM2017G105.DS		DM2017G160.DS
	3.0	DM3014G100.DS		DM3014G155.DS	
	4.0	DM4012G100.DS		DM4012G155.DS	
ATTACHMENT	80,000 rpm	MSB2-WX MSBF45135 MSBF4512135 MSBF4510C135 MSBF4512C135	MSB2-WX MSBF45135 MSBF4512135 MSBF4510C135 MSBF4512C135	MSB2-WX MSBF45155 MSBF45H155 MSBF4512155	MSB2-WX MSBF45155 MSBF45H155 MSBF4512155

Slim Bore: 1.19 mm



Tool Type	Diameter (mm)	Shaft Length 175 mm	Shaft Length 180 mm	Shaft Length 195 mm
	2.0		DM2017G180.DS	
	3.0	DM3014G175.DS		DM3014G195.DS
	4.0	DM4012G175.DS		DM4012G195.DS
ATTACHMENT	80,000 rpm	MSB2-WX MSBF45155 MSBF45H155 MSBF4512155	MSB2-WX MSBF45155 MSBF45H155 MSBF4512155	MSB2-WX MSBF45175 MSBF45H175 MSBF4512175 MSBF4512C175

Single-use Match-head Stainless-steel Slim Bur

Slim Bore: 1.19 mm



Tool Type	Diameter (mm)	Shaft Length 105 mm	Shaft Length 160 mm	Shaft Length 180 mm
	2.0	SM2002G105.DS	SM2002G160.DS	SM2002G180.DS
ATTACHMENT	80,000 rpm	MSB2-WX	MSB2-WX	MSB2-WX
		MSBF4580	MSBF45135	MSBF45155
		MSBF451280	MSBF4512135	MSBF45H155
			MSBF4510C135	MSBF4512155
			MSBF4512C135	

Telescoping Tube

Straight



Tool Type	Effective Length (mm)	Tube Diameter (mm)	Angle 0°	Angle 0° Protection Guard
	80	4.5	MSBF4580	
	135	4.5	MSBF45135	
	155	4.5	MSBF45155	
	155	4.5		MSBF45H155
	175	4.5	MSBF45175	
	175	4.5		MSBF45H175
ATTACHMENT	80,000 rpm	MSB2-WX	MSB2-WX	MSB2-WX

Curved



Tool Type	Effective Length (mm)	Tube Diameter (mm)	Angle 10°	Angle 12°
	135	4.5	MSBF4510C135	MSBF4512C135
	175	4.5		MSBF4512C175
ATTACHMENT	80,000 rpm	MSB2-WX	MSB2-WX	MSB2-WX

Angled



Tool Type	Effective Length (mm)	Tube Diameter (mm)	Angle 12°
	80	4.5	MSBF451280
	135	4.5	MSBF4512135
	155	4.5	MSBF4512155
	175	4.5	MSBF4512175
ATTACHMENT	80,000 rpm	MSB2-WX	

Single-use Spinal Bur

Round Diamond



Tool Type	Tube Diameter (mm)	Head Diameter (mm)	Effective Length 290 mm	Effective Length 350 mm
	3.0	3.0	JZTM29030J.DS	JZTM35030J.DS
	3.5	3.5	JZTM29035J.DS	JZTM35035J.DS
	4.0	4.0	JZTM29040J.DS	JZTM35040J.DS
ATTACHMENT		40,000 rpm	MSB1-ZJ	MSB1-ZJ

Round Stainless-steel



Tool Type	Tube Diameter (mm)	Head Diameter (mm)	Effective Length 290 mm	Effective Length 350 mm
	3.0	3.0	JZTM29030Q.DS	JZTM35030Q.DS
	3.5	3.5	JZTM29035Q.DS	JZTM35035Q.DS
	3.0	2.0	JZTM29030QH.DS	JZTM35030QH.DS
	3.5	2.5	JZTM29035QH.DS	JZTM35035QH.DS
ATTACHMENT		40,000 rpm	MSB1-ZJ	MSB1-ZJ

Cylinder Stainless-steel



Tool Type	Tube Diameter (mm)	Head Diameter (mm)	Effective Length 290 mm	Effective Length 350 mm
	3.0	2.0	JZTM29030QCH.DS	JZTM35030QCH.DS
	3.5	2.5	JZTM29035QCH.DS	JZTM35035QCH.DS
	3.0	2.0	JZTM29030QCH1.DS	JZTM35030QCH1.DS
	3.5	2.5	JZTM29035QCH1.DS	JZTM35035QCH1.DS
ATTACHMENT		40,000 rpm	MSB1-ZJ	MSB1-ZJ

Single-use Spinal Bur

Round Diamond



Tool Type	Tube Diameter (mm)	Head Diameter (mm)	Effective Length 290 mm	Effective Length 350 mm
	3.0	3.5	JZT29030J.DS	JZT35030J.DS
ATTACHMENT		25,000 rpm	JMSB1	JMSB1

Round Stainless-steel



Tool Type	Tube Diameter (mm)	Head Diameter (mm)	Effective Length 290 mm	Effective Length 350 mm
	3.0	3.5	JZT29030Q.DS	JZT35030Q.DS
ATTACHMENT		25,000 rpm	JMSB1	JMSB1

Cylinder Stainless-steel



Tool Type	Tube Diameter (mm)	Head Diameter (mm)	Effective Length 290 mm	Effective Length 350 mm
	3.0	2.0	JZT29030QCH1.DS	JZT35030QCH1.DS
ATTACHMENT		25,000 rpm	JMSB1	JMSB1

Single-use Bendable Spinal Bur

Round Diamond



2 sterile 5

Tool Type	Angle (°)	Tube Diameter (mm)	Head Diameter (mm)	Effective Length 290 mm	Effective Length 350 mm
	0~36	3.5	3.0	LB29030J.DS	LB35030J.DS
	0~36	3.5	3.5	LB29035J.DS	LB35035J.DS
		ATTACHMENT	25,000 rpm	JMSB1	JMSB1

Round Diamond - Coarse



2 sterile 5

Tool Type	Angle (°)	Tube Diameter (mm)	Head Diameter (mm)	Effective Length 290 mm	Effective Length 350 mm
	0~36	3.5	3.5	LB29035J1.DS	LB35035J1.DS
	0~36	3.5	5.0	LB29050J1.DS	
	0~36	3.5	6.0	LB29060J1.DS	
		ATTACHMENT	25,000 rpm	JMSB1	JMSB1

Single-use Bendable Spinal Bur

Round Stainless-steel



② sterile 5

Tool Type	Angle (°)	Tube Diameter (mm)	Head Diameter (mm)	Effective Length 290 mm	Effective Length 350 mm
	0~36	3.5	3.0	LB29030Q.DS	LB35030Q.DS
	0~36	3.5	3.5	LB29035Q.DS	LB35035Q.DS
		ATTACHMENT	25,000 rpm	JMSB1	JMSB1

Cylinder Stainless-steel



② sterile 5

Tool Type	Angle (°)	Tube Diameter (mm)	Head Diameter (mm)	Effective Length 290 mm	Effective Length 350 mm
	0~36	3.5	3.5	LB29035QC.DS	LB35035QC.DS
		ATTACHMENT	25,000 rpm	JMSB1	JMSB1

Single-use Nasal Shaver

Serrated Blade



Picture	Effective Length (mm)	Angle (°)	Tube Diameter 3.0 mm	Tube Diameter 3.6 mm	Tube Diameter 4.0 mm
	110	0	BPD11030A.DS	BPD11036A.DS	BPD11040A.DS
	110	12			BPD1104012A.DS
	110	40			BPD1104040A.DS
	110	60			BPD1104060A.DS
	110	90		BPD1103690A.DS	
	110	120		BPD11036120A.D S	
ATTACHMENT		6,000 rpm	PSB3	PSB3	PSB3
		6,000 rpm	PSB5	PSB5	PSB5

Serrated Blade



Picture	Effective Length (mm)	Angle (°)	Tube Diameter 4.0 mm
	110	0	BPD11040B.DS
ATTACHMENT		6,000 rpm	PSB3
		6,000 rpm	PSB5

Inferior Turbinate Blade



Picture	Effective Length (mm)	Angle (°)	Tube Diameter 3.0 mm
	110	0	BPD11030CH.DS
ATTACHMENT		6,000 rpm	PSB3
		6,000 rpm	PSB5

Single-use Nasal Bur

Round Diamond Bur



BZT1152515JH.DS

sterile

Picture	Effective Length (mm)	Angle (°)	Head Diameter 2.5 mm	Head Diameter 4.0 mm
	115	15	BZT1152515J.DS BZT1152515JH.DS	BZT1154015J.DS
	115	30	BZT1152530J.DS	BZT1154030J.DS
	120	70		BZT1204070J.DS
ATTACHMENT		13,000 rpm 13,000 rpm	PSB3 PSB5	PSB3 PSB5

Round Stainless-Steel Bur




sterile

Picture	Shaft Length (mm)	Angle (°)	Head Diameter 4.0 mm
	115	0	BZT11540QH.DS
	110	15	BZT1104015QH.DS
ATTACHMENT		13,000 rpm 13,000 rpm	PSB3 PSB5

Single-use Larynx Shaver

Serrated Blade



Picture	Effective Length (mm)	Angle (°)	Tube Diameter 4.0 mm
	225	18	HPD2254018A.DS
	275	18	HPD2754018A.DS

ATTACHMENT	6,000 rpm	PSB3	PSB3
	6,000 rpm	PSB5	PSB5

Single-use Adenoid & Tonsil Shaver

Serrated Blade



Picture	Effective Length (mm)	Angle (°)	Tube Diameter 4.0 mm
	110	12	YPD1104012A.DS
	110	40	YPD1104040A.DS

ATTACHMENT	6,000 rpm	PSB3	PSB3
	6,000 rpm	PSB5	PSB5

Single-use Arthroscopy Shaver

Serrated Blade



2 sterile 5

Picture	Effective Length (mm)	Tube Diameter 4.0 mm	Tube Diameter 4.5 mm	Tube Diameter 5.0 mm
	130	GPD13040A.DS	GPD13045A.DS	GPD13050A.DS
ATTACHMENT	6,000 rpm	PSB2	PSB2	PSB2

Serrated Blade



2 sterile 5

Picture	Effective Length (mm)	Tube Diameter 4.0 mm	Tube Diameter 4.5 mm	Tube Diameter 5.0 mm
	130	GPD13040B.DS	GPD13045B.DS	GPD13050B.DS
ATTACHMENT	6,000 rpm	PSB2	PSB2	PSB2

Serrated Blade



2 sterile 5

Picture	Effective Length (mm)	Tube Diameter 4.0 mm	Tube Diameter 4.5 mm	Tube Diameter 5.0 mm
	130	GPD13040C.DS	GPD13045C.DS	GPD13050C.DS
ATTACHMENT	6,000 rpm	PSB2	PSB2	PSB2

Double-edged Blade



2 sterile 5

Picture	Effective Length (mm)	Tube Diameter 4.0 mm	Tube Diameter 4.5 mm
	130	GPD13040B1.DS	GPD13045B1.DS
ATTACHMENT	6,000 rpm	PSB2	PSB2

Single-use Arthroscopy Bur

Round Stainless-Steel



② sterile 5

Picture	Effective Length (mm)	Tube Diameter 4.0 mm	Tube Diameter 5.0 mm	Tube Diameter 5.5 mm
	130	GZT13040Q.DS	GZT13050Q.DS	GZT13055Q.DS

ATTACHMENT 12,000 rpm PSB2

Round Stainless-Steel



② sterile 5

Picture	Effective Length (mm)	Tube Diameter 4.0 mm	Tube Diameter 5.0 mm	Tube Diameter 5.5 mm
	130	GZT13040QH.DS	GZT13050QH.DS	GZT13055QH.DS

ATTACHMENT 12,000 rpm PSB2 PSB2 PSB2

Cylinder Stainless-Steel



② sterile 5

Picture	Effective Length (mm)	Tube Diameter 4.0 mm	Tube Diameter 5.0 mm	Tube Diameter 5.5 mm
	130	GZT13040QC.DS	GZT13050QC.DS	GZT13055QC.DS

ATTACHMENT 12,000 rpm PSB2 PSB2 PSB2

Cylinder Stainless-Steel



② sterile 5

Picture	Effective Length (mm)	Tube Diameter 4.0 mm	Tube Diameter 5.0 mm	Tube Diameter 5.5 mm
	130	GZT13040QCH.DS	GZT13050QCH.DS	GZT13055QCH.DS

ATTACHMENT 12,000 rpm PSB2 PSB2 PSB2

Single-use Saw Blade

Small Bone - Saggital

② sterile 5

Picture	Effective Length (mm)	Tooth Width (mm)	Tooth Thickness 0.4 mm	Tooth Thickness 0.7 mm
	18	4	BJP180403.DS ● ●	
	31	10	BJP311003.DS ●	
	50	10		BJP501004.DS ●
ATTACHMENT		20,000 cpm	XPBJ ● XCBJ ●	XPBJ ●

Large Bone - Saggital

② sterile 5

Picture	Effective Length (mm)	Tooth Width (mm)	Thickness (mm)	Tooth Thickness 0.9 mm
	70	19	0.9	BJP701909.DS
	90	13	0.9	BJP901309.DS
	90	19	0.9	BJP901909.DS
ATTACHMENT		15,000 cpm	JJT1	

Large Bone-Sheathed

② sterile 5

Picture	Effective Length (mm)	Tooth Width (mm)	Thickness (mm)	Tooth Thickness 0.9 mm
	114	14	1.2	BJPH881412A.DS
	114	18	1.2	BJPH881812.DS
	130	18	1.2	BJPH1051812.DS
ATTACHMENT		15,000 cpm	JJT4	

Single-use Saw Blade

Small Bone - Reciprocating



Picture	Tooth Length (mm)	Tooth Thickness (mm)	Blade Thickness 0.6 mm
	15	1	WJP1506.DS
	24	1	WJP2406.DS
	31	1	WJP3106.DS
ATTACHMENT		20,000 cpm	XWFJ

Small Bone - Reciprocating



Picture	Length (mm)	Tooth Length (mm)	Blade Thickness 0.3 mm
	81	7	WJPZ81070304.DS
	85	23	WJPZ85230306.DS
	85	27	WJPZ85270306.DS
	91	33	WJPZ91330306.DS
	105	27	WJPZ105270306.DS
	111	33	WJPZ111330306.DS
ATTACHMENT		20,000 cpm	XWFJ

Large Bone - Reciprocating



Picture	Length (mm)	Tooth Thickness (mm)	Blade Thickness 0.6 mm
	60	1	WJP6006.DS
ATTACHMENT		15,000 cpm	JJT2

Single-use Saw Blade

Small Bone - Oscillating

2 **sterile** 5

Picture	Length (mm)	Tooth Width (mm)	Stop Depth (mm)	Angle (°)	Swing Distance 0.3 mm
	72	10	8	90	CBJP72100308.DS
	72	16	12	90	CBJP72160312.DS
	105	9	7	110	CBJP105090307110.DS
	107	9	12	110	CBJP107090312110.DS

ATTACHMENT 20,000 cpm XCBJ1

Single-use Twist Drill

Small Bone




② sterile 5

Tool Type	Head Diameter (mm)	Length 35 mm	Length 60 mm	Length 100 mm
	1.1	GT3511.DS		
	1.5	GT6015.DS		
	2.0	GT10020.DS		
	2.5	GT11025.DS		
ATTACHMENT	3,000 rpm	XKXZ2	XKXZ2	XKXZ2

Small Bone



② sterile 5

Tool Type	Head Diameter (mm)	Stop Depth (mm)	Length 130 mm
	1.1	4.0	GT130114B.DS
	1.1	6.0	GT130116B.DS
	1.5	4.0	GT130154B.DS
	1.5	6.0	GT130156B.DS
	1.5	8.0	GT130158B.DS
	1.8	6.0	GT130186B.DS
	1.8	8.0	GT130188B.DS
ATTACHMENT	80,000 rpm 40,000 rpm		MSB2A-W125 MSB1A-Z125 MSB-FJ1C

Single-use Twist Drill

Large Bone



Tool Type	Head Diameter (mm)	Length 145 mm	Length 155 mm	Length 195 mm
	3.2	GT14532.DS		
	4.0	GT15540.DS		
	4.5	GT19545.DS		
	6.0	GT19560.DS		
ATTACHMENT	1,500 rpm 1,500 rpm	ZJT14 ZJT15	ZJT14 ZJT15	ZJT14



CHONGQING XISHAN SCIENCE & TECHNOLOGY CO.,LTD.

Address: Jupiter S&T Development Center, No.9 Huangshan RD.(Middle),
High-tech Park, New North Zone, Chongqing

Tel:+86 - 23 - 68692230

Fax:+86 - 23 - 63211079

Website: english.xishantech.com.cn

Email: info@xishantech.com.cn